# DIASTAT<sup>®</sup> ANA

## (anti-nuclear antibody)

For professional use only  
Usage reserve aux professionnels  
Sólo para uso profesional  
Nur für den fachgebrauch  
Solo per uso professionale  
Apenas para utilização por profissionais  
Kun til professionel brug  
Endast för professionell användning

Document No. E-23-0103-08  
December, 2013

### DIASTAT<sup>®</sup> ANA

English:	page .....	2
Français:	page .....	13
Español:	página .....	24
Deutsch:	seite.....	35
Italiano:	pagina .....	46
Português:	pagina .....	57
Dansk:	side.....	68
Svenska:	sida.....	79

## **ENGLISH : INTENDED USE**

---

The DIASTAT<sup>®</sup> anti-nuclear antibody (ANA) test is a qualitative enzyme-linked immunosorbent assay (ELISA) for the detection of ANAs in human serum or EDTA plasma. It detects ANAs against Sm, Sm/RNP, Ro (SS-A), La (SS-B), Scl-70, Jo-1, dsDNA, histone and centromere antigens. The test may be used to screen out samples which are negative for all ANAs. Samples that give a positive test result should be further tested to identify the antigen-specific antibody or antibodies present. Quantitative/qualitative ELISA kits are available for the individual detection of Sm, Sm/RNP, Ro (SS-A), La (SS-B), dsDNA and centromere ANAs. Qualitative kits are available for the individual detection of Scl-70 and Jo-1. ANA detection represents one parameter in a multicriterion diagnostic process.

## **INTRODUCTION**

---

Systemic rheumatic diseases are autoimmune disorders such as systemic lupus erythematosus and lupus variants, polymyositis, Sjögren's syndrome, scleroderma, mixed connective tissue disease, CREST and rheumatoid arthritis. A general feature of systemic rheumatic diseases is the presence of circulating antibodies to a variety of cellular antigens. The detection and serological characterisation of specific autoantibodies plays an important role in the differential diagnosis of these diseases. Autoantibodies to nuclear antigens (ANAs) are a group of antibodies specific for nuclear antigens including Sm, Sm/RNP, Ro (SS-A), La (SS-B), Scl-70, Jo-1, dsDNA, histone and centromere antigens. A positive ANA test result provides presumptive evidence for systemic rheumatic disease<sup>1,2</sup>; further definition of specific antibody profiles is a valuable aid in the diagnostic process. A number of techniques are available for ANA detection including haemagglutination, immunodiffusion and indirect immunofluorescence<sup>3,4,5</sup>. These tests have inherent drawbacks - lack of sensitivity and reproducibility, and are unsuitable for large batch testing. ELISAs offer sensitivity, reproducibility, objectivity over other methods and flexibility in batch size testing. The DIASTAT<sup>®</sup> test detects total IgG and IgM antibodies to the extractable nuclear antigens Sm, Sm/RNP, Ro (SS-A), La (SS-B), Scl-70, Jo-1, dsDNA, histone and centromere antigens. It also detects serum autoantibodies producing homogenous, centromere and speckled/nucleolar patterns in conventional IFA testing. The DIASTAT<sup>®</sup> ELISA is an appropriate first line screen for total antibody. A negative result reduces, but does not exclude, the likelihood of systemic rheumatic disease. Positive results should be further analysed to identify the specific autoantibodies present.

## **PRINCIPLE OF THE ASSAY**

---

The wells of the microtitre strips are coated with highly purified HEp-2 cell extract. During the first incubation, specific ANAs in diluted serum or EDTA plasma bind to the antigen-coated surface. The wells are then washed to remove unbound components. In the second incubation, the Conjugate, enzyme-labelled antibodies to human IgG and IgM, binds any surface-bound autoantibodies. After further washing, specific autoantibodies are traced by incubation with the Substrate. Addition of Stop Solution terminates the reaction, resulting in a coloured end-product. The amount of Conjugate bound is measured in absorbance units and compared with that bound by the Reference Control.

## KIT COMPONENTS

<b>A</b>	IgG/IgM Conjugate	1 × 15mL	Alkaline phosphatase-labelled antibodies to human IgG and IgM, Tris buffer, protein stabiliser, <0.1% (w/v) sodium azide. <b>Ready-to-use.</b>	
<b>B</b>	Substrate	1 × 15mL	Mg <sup>2+</sup> , phenolphthalein monophosphate (PMP), buffer solution. <b>Ready-to-use.</b> Do not expose to light during storage..	
<b>C</b>	Stop Solution	1 × 15mL	Sodium hydroxide, EDTA, carbonate buffer (pH >10). <b>Ready-to-use..</b>	
<b>D</b>	Wash Buffer Concentrate (16X)	2 × 25mL	Borate buffer, 0.4% (w/v) sodium azide. <b>Dilute before use..</b>	
<b>E</b>	HEp-2-Coated Wells and Strip Holder	12 × 8 well microtitre strips	Coated with HEp-2 cell extract, in a resealable foil pack with desiccant.	
<b>F</b>	Sample Diluent 2 Concentrate (5X)	1 × 25mL	Phosphate buffer, protein stabiliser, 5% (w/v) Triton X-100, 0.5% (w/v) sodium azide. <b>Dilute before use.</b>	
<b>6</b>	ANA Reference Control	1 × 1.5mL	Human plasma, buffer, <0.1% (w/v) sodium azide. <b>Ready-to-use.</b>	
<b>+/-</b>	Positive Control Negative Control	1 × 0.2mL 1 × 0.1mL	Human plasma, <0.1% (w/v) sodium azide. <b>Dilute 1:101 with diluted Sample Diluent 2 before use, as for samples.</b>	
	Pack Leaflet			

## STORAGE OF REAGENTS

### *Opened Kit Stability*

A kit was opened and reused on three occasions over a three month period with no adverse effect on kit performance.

### *Handling and Procedural Notes*

1. Store kit components at 2-8°C and use until the expiry date on the labels. Do not use expired reagents.
2. Do not mix different lot numbers.
3. Do not freeze kits.
4. Wash Buffer Concentrate, Sample Diluent 2 Concentrate and Positive and Negative Controls must be diluted before use. All other reagents are ready-to-use.
5. Diluted Wash Buffer and diluted Sample Diluent 2 are stable at 2-8° C for up to 6 months if microbial contamination is avoided.
6. Replace surplus microtitre strips in the foil pack and store with the desiccant at 2-8° C, until required.
7. Do not expose Substrate to light during storage.
8. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.

### *Indications of Deterioration*

The Substrate should be pale yellow in colour. Pink colouring indicates contamination and the reagent must be discarded. Turbidity or precipitation in any component indicates deterioration and the component should be discarded.

### Sample Collection and Storage

The assay is recommended for serum/EDTA plasma samples; do not use lipaemic, haemolysed or turbid samples. Thoroughly mix thawed samples before assay and avoid repeated freeze/thawing. Do not heat-inactivate samples, this may yield false positive results.

Samples may be stored at or below -20° C, or assayed within 12 hours of collection.

## WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only.

### Safety Precautions

1. Adhere strictly to the instructions in this booklet, particularly for handling and storage conditions.
2. The Controls contain human plasma tested by FDA-cleared assays for hepatitis B surface antigen, HCV, HIV antigen and HIV antibodies and found to be non-reactive/negative. As no known test offers complete assurance that infectious agents are absent, Controls should be considered potentially infectious and handled with the same precautions as any other potentially biohazardous material. The CDC/NIH Health Manual "Biosafety in Microbiological and Biomedical Laboratories", 3rd edition, 1993, describes how these materials should be handled in accordance with Good Laboratory Practice. This is applicable in the USA.
3. Do not pipette by mouth.
4. Do not smoke, eat, drink or apply cosmetics in areas where kits and samples are handled.
5. Any skin complaints, cuts, abrasions and other skin lesions should be suitably protected.
6. The Controls, Conjugate, Sample Diluent 2 Concentrate and Wash Buffer Concentrate contain sodium azide which can react with lead and copper plumbing to form highly explosive metal azides. On disposal, drain with large quantities of water to prevent azide build-up.
7. The Stop Solution contains sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped up with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.
8. The substrate contains PMP, Bronidox L and Diethanolamine. Avoid contact with skin, eyes and respiratory system. If contact with skin, eyes or respiratory system occurs, rinse with water and seek medical advice.
9. Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

**Warning**

**B.**

**SUBS**

Contains: Diethanolamine

- |                 |  |
|-----------------|--|
| H319:           | Causes serious eye irritation.   |
| P264:           | Wash hands thoroughly after handling.  |
| P280:           | Wear protective gloves/protective clothing/eye protection/face protection.   |
| P305+P351+P338: | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313:      | If eye irritation persists: Get medical advice/attention.  |

C.

SOLN	STOP
------	------

**Warning**

Contains: Sodium hydroxide

- H315: Causes skin irritation.  
 H319: Causes serious eye irritation.  
 P264: Wash hands thoroughly after handling.  
 P280: Wear protective gloves/protective clothing/eye protection/face protection.  
 P302+P352: IF ON SKIN: Wash with plenty of soap and water.  
 P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P332+P313: If skin irritation occurs: Get medical advice/attention.  
 P337+P313: If eye irritation persists: Get medical advice/attention.

D. and F.

BUF	WASH	16X
DIL	SPE	5X

**Warning**

Contains: Sodium azide

- H302: Harmful if swallowed.  
 EUH032: Contact with acids liberates very toxic gas.  
 H412: Harmful to aquatic life with long lasting effects.  
 P264: Wash hands thoroughly after handling.  
 P280: Wear protective gloves/protective clothing/eye protection/face protection.  
 P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.  
 P273: Avoid release to the environment.

## PREPARATION

### **Materials/Equipment Required but not Provided**

1. 96 well plate/strip reader with 550nm filter (540-565nm is acceptable).
2. Precision pipettes to dispense 10µL, 100µL, 1mL. Automatic pipette to dispense 100µL. Automatic pipette to dispense 200µL for manual washing, automatic plate washer optional.
3. Glass/plastic measuring cylinders: 1×100mL, 1×400mL.
4. 1mL volume containers.
5. Distilled/deionised water.
6. Paper towels.
7. Timer for 30 and 60 minute intervals.

**Preparation for the Assay**

Allow all kit components, including the microtitre strips, to warm up to 18-25° C for 30-60 minutes before use. Mix reagents by gentle inversion.

**Do not dilute the Reference Control.**

Dilute the following reagents and mix thoroughly.

Reagent	Volume	Add
Wash Buffer Concentrate	1 vial	375mL distilled/deionised water
Sample Diluent 2 Concentrate	1 vial	100mL distilled/deionised water
Positive and Negative Controls/samples	10µL	1mL diluted Sample Diluent 2

Calculate the number of microtitre strips required for the current assay, and retain these in the microtitre strip holder. Return surplus strips to the foil pack with the desiccant, and store in the resealable plastic bag at 2-8°C until required. Ensure that all strips are securely held within the microtitre strip holder, with the assay identification tab lying along the bottom edge below row H.

**ASSAY PROTOCOL**

- Reference wells for identification.
- Pipette 100µL Reference Control in duplicate, pre-diluted Positive and Negative Controls, and pre-diluted patient samples into appropriate wells. Remember to change pipette tips between additions. This step should not exceed 15 minutes for any one set of Controls/samples.
- Incubate 60±10 minutes at 18-25° C.
- Decant strip contents by quick inversion over a sink suitable for the disposal of biological materials, bearing in mind the potential infective hazard of the samples. Blot inverted strips well with paper towels.
- Wash wells **three times** with a minimum of 200µL diluted Wash Buffer. **Decant and blot after each wash step.**
- Add 100µL IgG/IgM Conjugate to each well.
- Incubate 30±5 minutes at 18-25° C.
- Repeat steps 4 and 5.
- Add 100µL Substrate to each well.
- Incubate 30±5 minutes at 18-25° C. **Do not decant.**
- Add 100µL Stop Solution to each well, in the same order and rate as the Substrate. Tap wells gently to mix.
- Read strips within 24 hours at 550nm (540-565nm).

## CALCULATION AND INTERPRETATION OF RESULTS

**Consider each assay separately when calculating and interpreting results.**

Calculate the absorbance value (optical density) ratio for the Positive and Negative Controls, and for each sample.

$$\text{Absorbance Ratio} = \frac{\text{Sample or Control Absorbance Value}}{\text{mean Reference Control Absorbance Value}}$$

Users should calculate a cut-off between positive and negative samples that is specific to their patient populations. Results from the patient populations used in the Euro Diagnostica clinical trial suggest the following cut-off:

<u>Absorbance Ratio</u>	<u>Result Interpretation</u>
<0.7	<b>Negative</b>
≥0.7 to <1.0	<b>Borderline - recommend repeat testing on subsequent sample</b>
≥1.0	<b>Positive</b>

## QUALITY CONTROL

Ensure that adequate maintenance and calibration of the plate-reader is performed according to the manufacturer's instructions, and that the correct wavelength is employed.

Users should ensure they are fully acquainted with the instructions for the assay, particularly the Warnings and Precautions section, and the Handling and Procedural Notes. Users should demonstrate that they can obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results. It is recommended that the pre-diluted Positive and Negative Controls are run in duplicate in all assays to monitor the quality of the test procedure. Run the ready-to-use Reference Control in duplicate in all assays.

Assuming the precision specifications described by the manufacturer are met, failure of any Control to meet the Control ratio specifications below renders the assay invalid and patient results should not be reported. The operator may repeat the assay, having reviewed their procedure, or contact the distributor/manufacturer. If repeating the assay, prepare a fresh dilution of each Control and sample. Laboratories may wish to include in-house controls in each assay run. Store such control material at or below -20° C and avoid repeat freeze/thaw cycles. Preservatives such as sodium azide at 0.1% (w/v) will not affect sample results.

Levels of analytes identified in particular diseases are those established by the manufacturer for specific populations, and may not necessarily mirror the literature. Incidence levels, their relationship to specific diseases, reference ranges, and appropriate cut-off points should all be calculated for the specific populations serviced by users.

<b>Control Ratio Specifications</b>	
$\frac{\text{Positive Control Absorbance}}{\text{Reference Control Absorbance}}$	see Positive Control label
$\frac{\text{Negative Control Absorbance}}{\text{Reference Control Absorbance}}$	<0.7

## EXPECTED VALUES

169 serum samples from asymptomatic apparently healthy donors with no history of autoimmune or rheumatic disease were assayed for IgG and IgM autoantibodies to nuclear antigens. This population was assembled across four European centres; results were compared against conventional immunofluorescence using HEp-2 cells as substrate. Using a cut-off ratio of 0.7, the following agreement was found.

		DIASTAT	
		+ve	-ve
IFA	+ve	1	1
	-ve	10	158

Overall agreement = 94%

### *Comparison with Immunofluorescence*

363 serum samples from patients with connective tissue disease (CTD) were compared with immunofluorescence on samples from patients with connective tissue disease at a number of European centres. Positivity or negativity on IFA was assigned according to the trial centre criteria.

		DIASTAT	
		+ve	-ve
IFA	+ve	204	32
	-ve	8	119

Sensitivity = 86%

Specificity = 94%

Overall agreement = 89%

### *Comparison with Immunofluorescence Patterns*

IFA Pattern	n	DIASTAT ANA Positive	% Agreement
Centromere	35	34	97
Nucleolar/Speckled	26	25	96
Homogenous	35	32	91
Speckled	61	53	87
Nucleolar	47	36	77

### Comparison with Immunofluorescence by Individual Autoantibodies

The assay was compared with immunofluorescence, with individual autoantibody specificities. These were determined using the specific DIASTAT<sup>®</sup> ELISA kits.

Antibody Specificity	DIASTAT Positive (ANA)	Positive (IFA)	% Agreement
Ro (SS-A)	10	10	100
La (SS-B)	6	6	100
Sm	9	9	100
Sm/RNP	20	20	100
Scl-70	38	38	100
dsDNA	15	16	94
Centromere	61	62	98
Histone	12	13	92

28 samples that were identified as Jo-1 positive using the DIASTAT<sup>®</sup> anti-Jo-1 ELISA were run in the DIASTAT<sup>®</sup> ANA kit with the following results:

24/28 (85.7%) were positive in the ANA test.

3/28 (10.7%) fell within the borderline area. Their ratios were 0.99, 0.86, and 0.78.

1/28 (3.6%) gave a ratio of 0.67.

## PERFORMANCE DATA

### Imprecision

- Intra-assay imprecision** determined by testing three controls in three assays, using three operators and three kit batches, with replication of eight.

Control	Absorbance Value	SD	%CV
1	0.771	0.026	3.4
2	1.102	0.031	2.8
3	1.136	0.055	4.8

- Inter-assay imprecision** determined by testing six controls in 20 assays, using three operators and eight kit batches, with replication of four.

Control	Absorbance Value	SD	%CV
1	0.219	0.031	14.2
2	0.273	0.031	11.4
3	0.365	0.036	9.9
4	0.521	0.036	6.9
5	1.078	0.130	12.1
6	1.012	0.131	12.9

## LIMITATIONS OF USE

---

1. Although the presence of high titres of autoantibodies to nuclear antigens is indicative of systemic rheumatic diseases, the data must be considered in light of other clinical and laboratory findings.
2. Some individuals may have raised levels of ANAs with little or no evidence of clinical disease. By contrast, some patients with evidence of clinical disease may have undetectable levels of these antibodies.
3. The test does not identify the ANA antibody specificity. If a sample gives a positive result then further tests should be carried out to determine autoantibody specificity.
4. For repeat patient sampling, e.g. for monitoring, the same type of sample (serum or EDTA plasma) should be used throughout the study period.

## REFERENCES

---

1. Tan EM. Antinuclear Antibodies: Diagnostic Markers for Autoimmune Disease and Probes for Cell Biology. *Advances in Immunology*, **44**, 93-151, 1982.
2. Nakamura RM, Tan EM. Recent Advances in Laboratory Tests and the Significance of Autoantibodies to Nuclear Antigens in Systemic Rheumatic Diseases. *Clinics in Laboratory Medicine*, **6**, 41-53, 1986.
3. Sharp GC, Irvin WS, et al. Association of Autoantibodies to Different Nuclear Antigens with Clinical Patterns of Rheumatic Disease and Responsiveness to Therapy. *J Clin Invest*, **50**, 350-359, 1971.
4. Friou GJ. Antinuclear Antibodies: Diagnostic Significance and Methods. *Arthritis Rheum*, **10**, 151-159, 1967.
5. Eds. Nakamura RM, Greenwald CA, et al. Autoantibodies to Nuclear Antigens (ANA): Advances in Laboratory Tests and Significance in Systemic Rheumatic Diseases. American Society of Clinical Pathologists Press, Chicago, 1978.

## SUMMARY OF PROTOCOL

---

1. Dilute samples and Positive and Negative Controls 1:101. Do not dilute Reference Control.
2. Add 100µL of Reference Control in duplicate, pre-diluted Positive and Negative Controls and samples into referenced wells of the microtitre strip.
3. Incubate 60±10 minutes at 18-25° C.
4. Wash strips 3 times.
5. Add 100µL of IgG/IgM Conjugate to each well.
6. Incubate 30±5 minutes at 18-25° C.
7. Wash strips 3 times.
8. Add 100µL of Substrate to each well.
9. Incubate 30±5 minutes at 18-25° C.
10. Add 100µL of Stop Solution to each well.
11. Read absorbance at 550nm.

	Batch code / Code du lot / Código de lote / Chargenbezeichnung / Codice del lotto / Código do lote / Lotnummer / Satsnummer
	Catalogue number / Référence du catalogue / Número de catálogo / Bestellnummer / Numero di catalogo / Referência de catálogo / Bestillingsnummer / Katalognummer
	Use by / Utiliser jusque / Fecha de caducidad / Verwendbar bis / Utilizzare entro / Prazo de validade / Skal anvendes inden / Hållbar till
	Temperature limitation / Limites de temperature / Limite de temperatura / Zulässiger Temperaturbereich / Limite di temperatura / Limites de temperatura / Opbevares ved / Temperaturgrænser
	Biological risks / Risques biologiques / Riesgo biológico/Biogefährdung / /Rischio biologico / Risco biológico / Biologiske risici / Biologisk risk
	Consult instructions for use / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consultare le istruzioni per l'uso / Consulte as instruções de utilização / Læs brugsvejledning / Se bruksanvisning
	In Vitro Diagnostic Medical Device / Dispositif medical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / Dispositivo médico para diagnóstico <i>in vitro</i> / Medicinsk udstyr til in vitro-diagnostik / In vitro diagnostika
	Manufacturer / Fabricant / Fabricante / Hersteller / Fabbicante / Fabricante Producent / Tillverkare
	Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos/ Ausreichend für "n" Ansätze / Conteúdo suficiente para <n> ensaios / Contém o suficiente para "n" testes / Indeholder tilstrækkeligt til <n> tests / Innehållet tillräckligt för n tester
	Warning / Attention / Atención / Achtung / Attenzione / Atenção / Advarsel / Varning
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive / Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro / Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika / Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro / Conformidade com In Vitro Diagnostic Directive 98/79/EC, Diretiva de Dispositivos Médicos / Medicinsk udstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF / Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<b>CONJ</b>	Conjugate / Conjugué / Conjugado / Konjugat / Conjugato / Conjugado / Konjugat Konjugat
<b>SUBS</b>	Substrate / Substrat / Substrato / Substrat / Substrato / Substrato / Substrat / Substrat
<b>SOLN STOP</b>	Stop solution / Solution d'Arrêt / Solución de Parada / Stopplösung / Soluzione bloccante / Solução de paragem / Stopopløsning / Stopplösning
<b>BUF WASH 16 x</b>	Wash buffer concentrate (16 X) / Concentré tampon de (16X lavage)/ Concentrado de Búfer de lavado (16 X) / Washpuffer-Konzentrat (16 X) / Tampone di lavaggio concentrato (16 X) / Tampão de lavagem concentrado (16 X) / Vaskebufferkoncentrat (16 X) / Tvättbuffert koncentrat (16 X)
<b>Ag</b>	Hep-2-coated wells and strip holder / Cupules enduites de Hep-2 et Portebandes / Soporte para Bandas y Vasos Recubiertos con HEp-2 / Hep-2-beschichtete Vertiefungen und Streifenrahmen / Pozzetti rivestiti di Hep-2 e supporto per strip / Poços revestidos com HEp-2 e suporte de tiras / Hep-2-coatede brønde og strip- holder / Hep-2-klädda brunnar och striphållare
<b>DIL SPE 5 X</b>	Sample Diluent 2 Concentrate (5 X) / Concentré 2 diluent pour échantillons (5 X)/ Concentrado de Diluyente 2 de Muestra (5 X) / Probediluens 2-Konzentrat / Diluyente per campioni 2 concentrato (5 X) / Diluente de Amostras 2 concentrado (5 X) / Prøvediluent 2-koncentrat (5 X) / Provspädningsbuffert koncentrat (5 X)
<b>CONTROL REF</b>	ANA Reference Control / Témoin de référence ANA / Control de Referencia ANA / ANA Referenzkontrolle / Controllo di riferimento ANA / Controllo de referência de ANAs / ANA Referencekontrol / ANA referenskontroll
<b>CONTROL +</b>	Positive Controls / Témoins positives / Controles Positivos / Positiv-Kontrollen / Controlli Positivi / Controlos Positivos / Positive kontrollen / Positiva kontrollen
<b>CONTROL -</b>	Negative Controls / Témoins négatifs / Controles Negativos / Negativ-Kontrollen / Controlli negativi / Controlos Negativos / Negative kontrollen / Negativa kontrollen

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**EURO DIAGNOSTICA AB**